Package 'MALDIquantForeign'

August 11, 2014

```
Version 0.9
Date 2014-08-11
Title Import/Export routines for MALDIquant
Depends R (>= 3.0.0), methods, MALDIquant (>= 1.8)
Imports
     base64enc, digest, downloader, readBrukerFlexData (>= 1.7), readMzXmlData (>= 2.7), XML
Suggests knitr, testthat (>= 0.8), RNetCDF (>= 1.6.1)
Description This package reads (tab, csv, Bruker fid, Ciphergen XML,mzXML, mzML, imzML, Ana-
     lyze 7.5, CDF) and writes (tab, csv,msd, mzML) different file formats of mass spectrometry data
     into/from MALDIquant objects.
License GPL (>= 3)
URL http://strimmerlab.org/software/maldiquant/
     https://github.com/sgibb/MALDIquantForeign/
BugReports https://github.com/sgibb/MALDIquantForeign/issues/
LazyLoad yes
VignetteBuilder knitr
Author Sebastian Gibb [aut, cre], Pietro Franceschi [ctb]
Maintainer Sebastian Gibb <mail@sebastiangibb.de>
NeedsCompilation no
Repository CRAN
Date/Publication 2014-08-11 21:39:18
```

R topics documented:

MALDIquantForeign-package	2
export,AbstractMassObject-method	3
exportMsd,MassSpectrum-method	4
exportMzMl,MassSpectrum-method	5
exportTab,AbstractMassObject-method	6
import	8
importAnalyze	10
importBrukerFlex	11
importCdf	12
importCiphergenXml	13
importImzMl	14
importMzMl	15
importMzXml	16
importTxt	17
supportedFileFormats	18
	20

MALDIquantForeign-package

Import/Export routines for MALDIquant

Description

Index

This package reads and writes different file formats of mass spectrometry data into/from MALDIquant objects.

Details

Package: MALDIquantForeign

License: GPL (>= 3)

URL: http://strimmerlab.org/software/maldiquant/

Author(s)

Sebastian Gibb <mail@sebastiangibb.de>

References

http://strimmerlab.org/software/maldiquant/

```
export, AbstractMassObject-method Export\,files
```

Description

This function provides a general interface to export AbstractMassObject-class objects (e.g. MassSpectrum-class, MassPeaks-class) into different file formats.

Usage

```
## S4 method for signature 'AbstractMassObject'
export(x, file, type="auto", force=FALSE, ...)
## S4 method for signature 'list'
export(x, path, type, force=FALSE, ...)
```

Arguments

Х	a AbstractMassObject-class object or a list of AbstractMassObject-class objects.
file	character, file name.
path	character, path to directory in which the list of AbstractMassObject-class would be exported.
type	character, file format. If type is set to "auto" the file extension is used.
force	logical, If TRUE the file would be overwritten or path would be created.
	arguments to be passed to specific export functions.

Details

Specific export functions:

```
tab exportTab csv exportCsv msd exportMsd mzML exportMxMl
```

Author(s)

Sebastian Gibb

References

http://strimmerlab.org/software/maldiquant/

See Also

MassPeaks-class, MassSpectrum-class

Examples

exportMsd, MassSpectrum-method

Export to MSD files

Description

This function exports AbstractMassObject-class objects (e.g. MassSpectrum-class, MassPeaks-class) into mMass MSD files.

Usage

```
## S4 method for signature 'MassSpectrum'
exportMsd(x, file, force=FALSE, peaks, ...)
## S4 method for signature 'list'
exportMsd(x, path, force=FALSE, peaks, ...)
```

Arguments

```
x a MassSpectrum-class object or a list of MassSpectrum-class objects.

file character, file name.

path character, path to directory in which the list of AbstractMassObject-class would be exported.

peaks a MassPeaks-class object or a list of MassPeaks-class objects.

force logical, If TRUE the file would be overwritten or path would be created.

arguments to be passed to write.table.
```

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/,
mMass homepage: http://mmass.org/
```

See Also

MassPeaks-class, MassSpectrum-class

Examples

exportMzMl, MassSpectrum-method

Export to mzML files

Description

 $This function \ exports \ Abstract Mass Object-class \ objects \ (e.g. \ Mass Spectrum-class, Mass Peaks-class) into \ mzML \ files.$

Usage

```
## S4 method for signature 'MassSpectrum'
exportMzMl(x, file, force=FALSE, ...)
## S4 method for signature 'list'
exportMzMl(x, path, force=FALSE, ...)
```

Arguments

Х	a MassSpectrum-class object or a list of MassSpectrum-class objects.
file	character, file name.
path	character, path to directory in which the list of ${\tt AbstractMassObject-class}$ would be exported.
force	logical, If TRUE the file would be overwritten or path would be created.
	arguments to be passed to write.table.

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/, HUPO Proteomics Standards Inititative mzML 1.1.0 Specification: http://www.psidev.info/mzml_1_0_0
```

See Also

MassPeaks-class, MassSpectrum-class

Examples

```
exportTab, AbstractMassObject-method Export\ to\ text\ files
```

Description

This function exports AbstractMassObject-class objects (e.g. MassSpectrum-class, MassPeaks-class) into different text file formats.

Usage

```
## S4 method for signature 'AbstractMassObject'
exportTab(x, file, force=FALSE, ...)

## S4 method for signature 'list'
exportTab(x, path, force=FALSE, ...)

## S4 method for signature 'AbstractMassObject'
exportCsv(x, file, force=FALSE, ...)

## S4 method for signature 'list'
exportCsv(x, path, force=FALSE, ...)
```

Arguments

X	a AbstractMassObject-class object or a list of AbstractMassObject-class objects.
file	character, file name.
path	character, path to directory in which the list of AbstractMassObject-class would be exported.
force	logical, If TRUE the file would be overwritten or path would be created.
	arguments to be passed to write.table.

Details

```
exportTab and exportCsv use write.table with different defaults (sep="\t" in exportTab and sep="," in exportCsv).
```

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/
```

See Also

```
MassPeaks-class, MassSpectrum-class, write.table
```

8 import

```
## export a single spectrum
exportTab(s[[1]], file="spectrum.tab")

## export a list of spectra and use ; as separator
exportCsv(s, path="spectra", sep=";", force=TRUE)

## End(Not run)
```

import

Import files

Description

This function provides a general interface to import different file formats into MassSpectrum-class or MassPeaks-class objects.

Usage

```
import(path, type = "auto", pattern, excludePattern = NULL,
  removeEmptySpectra = TRUE, centroided = FALSE, massRange = c(0, Inf),
  minIntensity = 0, verbose = TRUE, ...)
```

Arguments

path character, path to directory or file which should be read in.

type character, file format. If type is set to "auto" MALDIquant tries to detect the

correct file type automatically. It often depends on the file extension (if path is a directory the most represented file extension is used; pattern argument is

ignored).

pattern character, a regular expression to find files in a directory (see details).

excludePattern character, a regular expression to exclude files in a directory (see details).

removeEmptySpectra

logical, should empty spectra excluded?

centroided logical, if centroided=FALSE (default) the data are treated as not centroided

and a list of MassSpectrum-class objects is returned. Use centroided=TRUE

to assume centroided data and get a list of MassPeaks-class objects.

massRange double, limits of mass import (left/minimal mass, right/maximal mass).

minIntensity double, minimal intensity to import.

verbose logical, verbose output?

... arguments to be passed to specific import functions.

import 9

Details

Specific import functions:

```
txt importTxt
tab importTab
csv importCsv
fid importBrukerFlex
ciphergen importCiphergenXml
```

mzXML importMzXml
mzML importMzMl
imzML importImzMl
analyze importAnalyze
cdf importCdf

path: In addition to the above mentioned file types the following (compressed) archives are supported, too: zip, tar, tar.gz, tar.bz2, tar.xz. The archives are uncompressed in a temporary directory. Afterwards the import function is called (with type="auto").

pattern: Sometimes unusual file extensions are used (e.g. "*.xml" for mzXML files). In this case a specific pattern could be defined to import files with an unusual file extension (e.g. pattern="^.*\.xml\$" to read all *.xml files in a directory; see regexp for details).

excludePattern: Sometimes some files should be excluded. E.g. to ignore additional aquired Bruker LIFT spectra (MALDI-TOF/TOF; which are not support, yet) you could use excludePattern="([[:digit:]\.]+)L

Value

a list of MassSpectrum-class or MassPeaks-class objects (depending on the centroided argument).

Author(s)

Sebastian Gibb

References

http://strimmerlab.org/software/maldiquant/

See Also

MassSpectrum-class, MassPeaks-class

```
library("MALDIquant")
library("MALDIquantForeign")

## get example directory
exampleDirectory <- system.file("exampledata", package="MALDIquantForeign")</pre>
```

10 importAnalyze

```
## import mzXML files
s <- import(exampleDirectory, type="mzXML")

## import tab delimited file with different file extension (default: *.tab)
s <- import(exampleDirectory, type="tab", pattern="^.*\\.txt")

## import single mzML file
s <- import(file.path(exampleDirectory, "tiny1.mzML1.1.mzML"))

## import gzipped csv file
s <- import(file.path(exampleDirectory, "compressed", "csv1.csv.gz"))</pre>
```

importAnalyze

Import Analyze 7.5 files

Description

This function imports files in Analyze 7.5 file format into MassSpectrum-class or MassPeaks-class objects.

Usage

```
importAnalyze(path, ...)
```

Arguments

path character, path to directory or file which should be read in.
... arguments to be passed to import.

Value

a list of MassSpectrum-class or MassPeaks-class objects (depending on the centroided argument).

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/
http://www.grahamwideman.com/gw/brain/analyze/formatdoc.htm, http://eeg.sourceforge.
net/ANALYZE75.pdf
```

See Also

MassSpectrum-class, MassPeaks-class

importBrukerFlex 11

importBrukerFlex

Import Bruker Daltonics *flex files

Description

This function imports files in Bruker Daltonics *flex-series file format into MassSpectrum-class or MassPeaks-class objects.

Usage

```
importBrukerFlex(path, ...)
```

Arguments

path character, path to directory or file which should be read in.
... arguments to be passed to readBrukerFlexFile.

Value

a list of MassSpectrum-class or MassPeaks-class objects (depending on the centroided argument).

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/
```

See Also

MassSpectrum-class, MassPeaks-class, readBrukerFlexFile

```
library("MALDIquant")
library("MALDIquantForeign")

## get example directory
exampleDirectory <- system.file("exampledata", package="MALDIquantForeign")
s <- importBrukerFlex(exampleDirectory)</pre>
```

12 importCdf

importCdf

Import CDF files

Description

This function imports files in NetCDF file format into MassSpectrum-class or MassPeaks-class objects.

Please note that the *RNetCDF* is needed.

Usage

```
importCdf(path, ...)
```

Arguments

```
path character, path to directory or file which should be read in.
... arguments to be passed to import.
```

Value

a list of MassSpectrum-class or MassPeaks-class objects (depending on the centroided argument).

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/
```

See Also

MassSpectrum-class, MassPeaks-class

```
library("MALDIquant")
library("MALDIquantForeign")

## get example directory
exampleDirectory <- system.file("exampledata", package="MALDIquantForeign")

## import
if (require("RNetCDF")) {
    s <- importCdf(exampleDirectory)
} else {
    message("You have to install the RNetCDF package to use importCdf.")
}</pre>
```

importCiphergenXml 13

importCiphergenXml

Import Ciphergen XML files

Description

This function imports files in Ciphergen XML file format into MassSpectrum-class or MassPeaks-class objects.

Usage

```
importCiphergenXml(path, ...)
```

Arguments

path character, path to directory or file which should be read in.
... arguments to be passed to import.

Value

a list of MassSpectrum-class or MassPeaks-class objects (depending on the centroided argument).

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/
```

See Also

MassSpectrum-class, MassPeaks-class

```
library("MALDIquant")
library("MALDIquantForeign")

## get example directory
exampleDirectory <- system.file("exampledata", package="MALDIquantForeign")

## import
s <- importCiphergenXml(exampleDirectory)</pre>
```

14 importImzMl

importImzMl

Import imzML files

Description

This function imports files in imzML file format into MassSpectrum-class or MassPeaks-class objects.

Usage

```
importImzMl(path, ...)
```

Arguments

```
path character, path to directory or file which should be read in.
... arguments to be passed to import.
```

Value

a list of MassSpectrum-class or MassPeaks-class objects (depending on the centroided argument).

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/,
Definition of imzML format: http://www.imzml.org/
```

See Also

MassSpectrum-class, MassPeaks-class

```
library("MALDIquant")
library("MALDIquantForeign")

## get example directory
exampleDirectory <- system.file("exampledata", package="MALDIquantForeign")

## import
s <- importImzMl(exampleDirectory)</pre>
```

importMzMl 15

importMzMl

Import mzML files

Description

This function imports files in mzML file format into MassSpectrum-class or MassPeaks-class objects.

Usage

```
importMzMl(path, ...)
```

Arguments

```
path character, path to directory or file which should be read in.
... arguments to be passed to import.
```

Value

a list of MassSpectrum-class or MassPeaks-class objects (depending on the centroided argument).

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/,
Definition of mzML format: http://www.psidev.info/mzml_1_0_0%20
```

See Also

MassSpectrum-class, MassPeaks-class

```
library("MALDIquant")
library("MALDIquantForeign")

## get example directory
exampleDirectory <- system.file("exampledata", package="MALDIquantForeign")

## import
s <- importMzMl(exampleDirectory)</pre>
```

16 importMzXml

importMzXml

Import mzXML files

Description

This function imports files in mzXML file format into MassSpectrum-class or MassPeaks-class objects.

Usage

```
importMzXml(path, ...)
```

Arguments

```
path character, path to directory or file which should be read in.
... arguments to be passed to readMzXmlFile.
```

Value

a list of MassSpectrum-class or MassPeaks-class objects (depending on the centroided argument).

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/,
Definition of mzXML format: http://tools.proteomecenter.org/mzXMLschema.php
```

See Also

```
MassSpectrum-class, MassPeaks-class, readMzXmlFile
```

```
library("MALDIquant")
library("MALDIquantForeign")

## get example directory
exampleDirectory <- system.file("exampledata", package="MALDIquantForeign")

## import
s <- importMzXml(exampleDirectory)</pre>
```

importTxt 17

importTxt

Import text files

Description

This function imports different text file formats into MassSpectrum-class or MassPeaks-class objects.

Usage

```
importTxt(path, ...)
importTab(path, ...)
importCsv(path, ...)
```

Arguments

```
path character, path to directory or file which should be read in.
... arguments to be passed to read.table.
```

Details

importTab, importTxt and importCsv use read.table with different defaults.

Value

a list of MassSpectrum-class or MassPeaks-class objects (depending on the centroided argument).

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/
```

See Also

MassSpectrum-class, MassPeaks-class, read.table

supportedFileFormats

Examples

```
library("MALDIquant")
library("MALDIquantForeign")

## get example directory
exampleDirectory <- system.file("exampledata", package="MALDIquantForeign")

## import txt files
s <- importTxt(exampleDirectory)

## import csv files
s <- importCsv(exampleDirectory)</pre>
```

supportedFileFormats Supported file formats

Description

This function prints all file formats supported by MALDIquantForeign-package.

txt

Usage

```
supportedFileFormats()
```

Details

Import:

tab	importTab
csv	importCsv
fid	importBrukerFlex
ciphergen	importCiphergenXml
mzXML	importMzXml
mzML	importMzMl
imzML	<pre>importImzMl</pre>

importTxt

imzML importImzMl analyze importAnalyze cdf importCdf

Export:

```
tab exportTab csv exportCsv msd exportMsd mzML exportMxMl
```

supportedFileFormats 19

Value

a list with two named elements (import and export) containing a character vector of supported file types.

Author(s)

Sebastian Gibb

References

```
http://strimmerlab.org/software/maldiquant/
```

See Also

```
export, import
```

```
library("MALDIquantForeign")
supportedFileFormats()
```

Index

```
*Topic package
    MALDIquantForeign-package, 2
                                                exportTab,AbstractMassObject-method,6
                                                exportTab,list-method
export, 19
                                                         (exportTab, AbstractMassObject-method),
export
        (export,AbstractMassObject-method),
                                                import, 8, 9, 10, 12-15, 19
export,AbstractMassObject-method, 3
                                                importAnalyze, 9, 10, 18
                                                importBrukerFlex, 9, 11, 18
export,list-method
        (export,AbstractMassObject-method),
                                                importCdf, 9, 12, 18
                                                importCiphergenXml, 9, 13, 18
                                                importCsv, 9, 18
exportCsv, 3, 18
                                                importCsv (importTxt), 17
exportCsv
        (exportTab, AbstractMassObject-method), importImzMl, 9, 14, 18
                                                importMzMl, 9, 15, 18
exportCsv,AbstractMassObject-method
                                                importMzXml, 9, 16, 18
        (exportTab, AbstractMassObject-method), importTab, 9, 18
                                                importTab(importTxt), 17
                                                importTxt, 9, 17, 18
exportCsv,list-method
        (exportTab,AbstractMassObject-method),
                                                MALDIquantForeign-package, 2
        6
exportMsd, 3, 18
                                                read.table.17
exportMsd
                                                readBrukerFlexFile. 11
        (exportMsd, MassSpectrum-method),
                                                readMzXmlFile, 16
                                                regexp, 9
exportMsd,list-method
        (exportMsd, MassSpectrum-method),
                                                supportedFileFormats, 18
exportMsd, MassSpectrum-method, 4
                                                write.table, 4, 6, 7
exportMzMl, 3, 18
exportMzMl
        (exportMzMl, MassSpectrum-method),
exportMzMl,list-method
        (exportMzMl, MassSpectrum-method),
exportMzMl, MassSpectrum-method, 5
exportTab, 3, 18
exportTab
        (exportTab,AbstractMassObject-method),
```